



TrioDocs


Version: 0.2.1

Date: May 08, 2025

Download the latest version on:  
<https://docs.diy-trio.org>

# Shortcuts

# Shortcuts



## Allow Bolusing with Shortcuts

**Default:** OFF

Enabling this setting allows the iOS Shortcuts App to send bolus commands to Trio.



### Tip

Disabling this setting will still allow other Shortcut commands, like Temp Targets, Add Carbs, and Start/End Overrides.

## Available Shortcut actions:

- [Bolus](#) (if enabled)
- [List State](#)
- [Add Carbs](#)
- [Activate Override](#)
- [Cancel Override](#)
- [Temporary Target](#)
- [Cancel Temporary Target](#)

- [Restart Live Activity](#)
- 

## Bolus

If [Allow Bolusing with Shortcuts](#) is enabled, you will be able to use shortcuts to deliver a bolus.

---

## List State

Shows the last glucose reading, glucose trend arrow, time since last reading, and glucose delta. It also shows current IOB and COB.

---

## Add Carbs

Adds a carb entry into Trio.



### Tip

You can only add a carb entry using shortcuts. You cannot add a fat or protein entry.

---

## Activate Override

Activates a chosen preset override.



### Tip

You must first save the Preset Override on the [Adjustments](#) screen before you can activate an override using Shortcuts.

---

## Cancel Override

Cancels an already set override.

---

## Temporary Target

Activates a chosen preset temp target.



### Tip

You must first save the Preset Temp Target on the [Adjustments](#) screen before you can activate a temp target using Shortcuts.

---

## Cancel Temporary Target

Cancels and already active temp target.

---

## Restart Live Activity

Restarts Trio's live activity.

---